



Glucose Prediction from Laser Speckle Image Patterns

~ Sailee Prashant Allyadwar

1. Project Overview

Objective

Develop a machine learning model that predicts blood glucose levels (mg/dL) from images of biological samples (e.g., test strips, skin speckle patterns).

Approach

- Extract texture features (GLCM, LBP) from grayscale images.
- Train regression models (XGBoost, Random Forest, etc.) to predict glucose levels.
- Optimize hyperparameters for best performance.

Key Challenges

- Limited dataset size (~300 images).
- High prediction errors (MAE: 65–87 mg/dL).
- Weak correlation (R^2 : 0.14–0.49).

2. Dataset

Structure

text

Copy

Download

F:\8.Glucometer Setup IITH\
├── train_data/
│ ├── 100 mg/dL/ # Folder per concentration
│ │ ├── image1.jpg
│ │ └── ...
│ ├── 200 mg/dL/
│ │ └── ...
└── test_data/
 └── ...

- Classes: Glucose concentrations (e.g., 100, 150, 200 mg/dL).
- Images: Grayscale, resized to 256×256 px.

Preprocessing

1. Validation
 - Check for corrupt/blank images.
 - Remove uniform-color images.
2. Normalization
 - Resize to 256×256.
 - Contrast normalization (cv2 . NORM_MINMAX).

3. Feature Extraction

Methods

Feature Type	Description	Parameters
GLCM	Gray-level co-occurrence matrix	Distances: [1, 3, 5], Angles: [0°, 45°, 90°, 135°]
LBP	Local Binary Patterns	Radius: [1, 3, 5], Points: 8×radius
Speckle Statistics	Mean, median, std, skewness, kurtosis	--
Histogram	Intensity distribution (16 bins)	Range: [0, 256]

1. Gray-Level Co-Occurrence Matrix (GLCM) Features

GLCM analyzes texture by calculating how often pairs of pixel intensities occur in an image at a given distance and angle.

Key Parameters

- Distance (d): 1, 3, 5 pixels
Defines how far apart pixel pairs should be.
- Angles (θ): 0°, 45°, 90°, 135°
Directions in which pixel pairs are analyzed.
- Properties Calculated:
 - Contrast: Measures local intensity variations.
 - High contrast = Sharp edges/textures.

$$\sum_{i,j} (i - j)^2 \cdot P(i, j)$$

- Formula:

- Dissimilarity: Similar to contrast but linear (less sensitive to outliers).

$$\sum_{i,j} |i - j| \cdot P(i, j)$$

- Formula:

- Homogeneity: Measures closeness of GLCM values to the diagonal.

- High homogeneity = Smooth textures.

$$\sum_{i,j} \frac{P(i, j)}{1 + |i - j|}$$

- Formula:

- Energy (Angular Second Moment): Measures uniformity of pixel pairs.

- High energy = Very similar pixels.

$$\sum_{i,j} P(i, j)^2$$

- Formula:

- Correlation: Measures linear dependency between pixels.

$$\sum_{i,j} \frac{(i - \mu_i)(j - \mu_j)P(i, j)}{\sigma_i \sigma_j}$$

- Formula:


Why GLCM?

- Captures texture patterns (e.g., smooth vs. grainy regions).
- Effective for biological samples where glucose affects speckle patterns.

Original:



Output:



2. Local Binary Patterns (LBP)

LBP encodes local texture by comparing each pixel to its neighbors.

Key Parameters

- Radius (R): 1, 3, 5 pixels
Defines neighborhood size.
- Points (P): $8 \times \text{radius}$
Number of neighboring pixels to compare.
- Method: 'uniform'
Reduces feature dimensions by grouping similar patterns.

How It Works

$$\text{LBP} = \sum_{p=0}^{P-1} s(g_p - g_c) \cdot 2^p$$

1. For each pixel, compare intensity with neighbors:
2. where $s(x) = 1$ if $x \geq 0$, else 0.
3. Compute histogram of LBP codes (bins = 0 to $P+2$).
4. Normalize histogram to sum to 1.

Why LBP?

- Rotation-invariant (uniform patterns).
- Captures micro-textures (e.g., speckle noise).

Original :



Output:



3. Speckle Statistics

Measures intensity distribution properties.

Features Extracted

1. Speckle Contrast


$$\frac{\sigma}{\mu} \quad (\text{Standard deviation} / \text{Mean})$$

- High contrast = More noise (e.g., higher glucose variability).
2. Intensity Mean/Median
 - Average/median pixel value.
3. Standard Deviation
 - Measures intensity spread.

4. Skewness
 - Asymmetry of intensity distribution.
 - Positive skew = More dark pixels.
5. Kurtosis
 - "Peakedness" of intensity distribution.
 - High kurtosis = Sharp intensity peaks.

Why Speckle Stats?

- Glucose changes may alter light scattering properties.



4. Models & Training

Algorithms Tested

Model	Tuning Method	Best MAE (mg/dL)	R²
XGBoost	Random Search	65.88 (Train)	0.4896
Random Forest	Grid Search	38.67 mg/dL(Train) 87.96 mg/dL(Validation) 82.20 mg/dL(Testing)	0.8527(Testing) -0.0301(Validation) 0.2255(Testing)

Random Forest Output:

Sample Predictions:

Predicted: 183 mg/dL		Actual: 150 mg/dL
Predicted: 237 mg/dL		Actual: 400 mg/dL
Predicted: 187 mg/dL		Actual: 300 mg/dL
Predicted: 283 mg/dL		Actual: 400 mg/dL
Predicted: 170 mg/dL		Actual: 250 mg/dL

XGBoost Output:

Sample Predictions:

Predicted: 184 mg/dL		Actual: 300 mg/dL
Predicted: 210 mg/dL		Actual: 250 mg/dL
Predicted: 201 mg/dL		Actual: 250 mg/dL
Predicted: 197 mg/dL		Actual: 200 mg/dL
Predicted: 181 mg/dL		Actual: 250 mg/dL

Hyperparameter Tuning (XGBoost)

```
python
Download
Best Parameters:
- learning_rate: 0.014
- max_depth: 3
```

- n_estimators: 109
- subsample: 0.825
- colsample_bytree: 0.872

Evaluation Metrics

Split	MAE (mg/dL)	RMSE (mg/dL)	R²
Training	65.88	79.32	0.4896
Validation	87.31	109.13	0.1388
Test	84.08	102.58	0.1725

5. Performance Analysis

Key Issues

1. High Error Rates
 - o MAE > 65 mg/dL is clinically unacceptable (target: <15 mg/dL).
 - o Poor performance on extremes (e.g., predicted 275 vs. actual 100 mg/dL).
2. Low R²
 - o Model explains <50% of variance.
3. Overfitting
 - o Training MAE (65.88) << Validation MAE (87.31).

Sample Predictions

Actual (mg/dL)	Predicted (mg/dL)	Error (%)
150	193	+28.7%
100	275	+175%
250	182	-27.2%

6. Recommendations for Improvement

Data-Level

1. Increase Dataset Size
 - o Collect 1,000+ images across glucose ranges.
2. Augmentation
 - o Rotate, flip, adjust contrast/blur to synthetically expand data.
3. Class Balancing
 - o Ensure equal samples per glucose level.

Model-Level

1. Try Simpler Models
 - o Linear Regression (baseline) or SVR with RBF kernel.
2. Adjust XGBoost
 - o Increase learning_rate (0.05–0.2).
 - o Use early stopping.
3. Feature Engineering
 - o Add color features (if RGB images available).
 - o Experiment with CNN-based feature extraction.

Evaluation

- 1. Clinical Accuracy
 - Report % predictions within $\pm 15\%$ of actual values.
- 2. Error Analysis
 - Visualize residuals vs. glucose levels.

7. Code Structure

Key Files

File	Purpose
<code>python6.py</code>	Main training/prediction script
<code>glucose_model.pkl</code>	Saved XGBoost model
<code>glucose_scaler.pkl</code>	Feature scaler

Usage

```
bash
Copy
Download
# Train model
python python6.py

# Predict single image
python predict.py --image test_image.png
```

8. Future Work

- Integrate deep learning (CNNs).
- Deploy as a mobile/web app.
- Validate with clinical trials.